

## InSync Universal Paste Stain & Glaze Quick Start

### Intended Use:

A Type 1 ceramic for the staining and glazing of the following restoration types: Zirconia, Titanium, Lithium Disilicate and Metal Ceramic.

### Pre-glaze Preparation:

Lightly air-abrade using aluminum oxide. Steam clean or place in ultrasonic prior to glazing.

### Before Use:

- ☑ **IMPORTANT:** Do not pour out excess liquid in glaze jar.
- ☑ **MIX THE STAIN AND GLAZE PASTE THOROUGHLY BEFORE USE EVERY DAY** to maintain consistency and proper handling properties throughout the life of the jar.
- ☑ **DO NOT** bring glaze into contact with water! Glaze material that has been contaminated by water will appear cloudy white after firing. Be sure to always use a clean, dry brush.

### Storage:

Store under dry conditions, +12°C to +38°C

### Firing Parameters - Stain & Glaze:

Fired to:	Entry Time (min)	Low Temp °C	Vac Start °C	Heat Rate °C/min	High Temp °C	Hold Time (min)
Full Contour Zirconia	6	450	620	45	760	1(air)
InSync Layering	6	450	600	45	730	1(air)
Press to Zirconia	6	450	600	45	760	1(air)
Lithium Disilicate	6	450	600	45	730	1(air)
Metal Ceramic	6	450	600	45	730	1(air)

CTE (25-500°C) [ $\times 10^{-6} \text{K}^{-1}$ ]  $\pm 0.5$ : Stain 7.5 (2x) | Glaze 7.3 (1x)

For long-span ZrO2 bridges, an increase in the firing temperature of 10°C to 20°C is recommended.

**Quickstart:** <http://bit.ly/InSyncQuickStart>



Pre-Glaze Preparation: Important pre-glaze preparation  
 One-Step Stain & Glaze: Stain & glaze in one firing  
 Two-Step Stain & Glaze: Stain & glaze in two separate firings