



Basic Powder Guidance

Vita Lumin shade guide

	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaque	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
OD	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Dentin	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Enamel	58	58	59	59	60	57	58,59	59	59	60	59	59	60	60	59	59
Neck Translucent*	52	52	53	53	53,55	52	52	52,54	53,54	51,52	52,53	53,54	52,55	52	53,54	52,54
Make In* (mix 1:1 w/Dentin)	61	62	63	63	65	62	62	63	65	61	64	64	66	64	66	66
SI*	01	02	03	04	04	01	02	03	04	02	02	03	05	03	03	04
Clears	CL-0, NT, & TI-4 are transparent powders that work well with most shades															

* HT, MI, SI, and TI powders are not shade specific - the materials indicated are merely suggestions of complimentary colors for the indicated shades. See "How to Apply the Unique Powders of the Creation Porcelain" system for more information.

Shoulder Porcelains

Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaque Bead	27	28	28	29	29	27	28	28	28	27	28	28	28	28	28	28
SP Blend	22;21	22	23;25	24;25	24;25	21;22	21;22	23	24;23	22;26	26;22	26;25	26;25	22;25	22;25;26	22;25
SP Blend Ratio	3:1	100%	2:1	1:1	2:1	4:1	2:1	100%	3:1	1:1	2:1	3:1	4:1	4:1	3:2:1	3:1

