

Imagine[®]

ZIRCONIA



Because beauty alone is not enough.

Strong. Reliable. Proven.

Now Available in Multilayered and Unshaded

Beautiful.



Imagine® mimics the translucency of natural teeth.

- More translucent than traditional Zirconia and as translucent as pressed ceramics
- Vital and accurate shading through customized dying liquid application with unshaded Imagine
- Seamless color gradation and out of the box accuracy with new 7 color multilayered Imagine zirconia
- Mimics a level of intraoral value like that of natural teeth, not the stark, higher values seen with traditional zirconia

Strong.

Imagine® is the strongest of all comparable Zirconia materials and 2X the strength of lithium disilicate.

- Exceeds minimum flexural strength for a class 4 Type II material up to 35% (ISO 6872:2015)
- Strong enough to use in the anterior and posterior region of the mouth & for bridges up to 3 units in length
- Can be prepped as thin as 0.8; enables dentist to preserve up to 30% of tooth structure compared to lithium disilicate

Reliable.

Imagine® consistently delivers the results your doctors demand.

- Careful selection of the finest raw materials and state of the art processing result in fewer remakes
- Highly controlled pressing and sintering techniques that guarantee disc to disc consistency

Proven.

Imagine® has a long track record...that's peace of mind for you and your dentist customers.

- Jensen was one of the first to receive 510k clearance to market this type of material
- Hundreds of thousands of restorations produced over 2 years of sales

Imagine Unshaded Zirconia available in a range of 12-22mm x 98mm discs & 14-18mm Lava™ frames.

Imagine Multilayered Zirconia available in 14 & 18mm x 98mm discs in 8 shades A1-A4, B1-B3.