



Preciso

Power

## NEW Identica T500 Scanner

Add more power, control and greater flexibility to your digital workflow with Jensen's new Identica T500 Scanner from Medit.

### Why choose the Identica T500?

#### POWERFUL

- Save time by scanning a full arch in 12 seconds or 8 single dies in 19 seconds
- Start scanning with the touch of a sensor
- Scan up to 16 single dies simultaneously
- Import/export Stl files anytime, anywhere
- Preciso CAD software powered by EXOCAD

#### ACCURATE


- Dual 2.0 MP cameras improve scan quality and details
- New impression scanning rivals scanned model in accuracy
- <7 micron scan accuracy

#### INNOVATIVE

- Process models with less scans and clicks; what used to take 7 steps now only takes 3
- 3 Axis impression arm allows for double sided impression scanning in one step
- Versatile grip plate holds what you need where you need it
- Full size articulator scanning
- Patent pending implant position scanning from an impression, no need to pour a plaster model
- New collab 2017 scan software offers never been seen features and interfaces

#### COST-EFFECTIVE

- Accessories included
- Auto start/stop technology extends scanner and bulb life while reducing electricity consumption
- Top tier performance and features, entry level pricing starts at \$14,495





Contact your Jensen Sales Consultant at 800.243.2000 to learn more about how Preciso can help improve quality and lower overall costs.



## ADVANCED DATA PROCESSING & HIGHER RESOLUTION CAMERAS

Identica T500 dramatically improves the most important aspects of scanner performance. It has revolutionized detail of data with the state-of-the-art software technology, higher resolution cameras and blue light scanning technology. The latest software technology provides you with the cleanest and sharpest data. It introduces an advanced new camera with dual 2.0 MP resolution offering the best performance and scan data quality ever in Identica scanners. For the ultimate in scan data quality, our new T500 presents brilliance of clarity in in your work and will save you precious time and money by minimizing trial and error to adjust restoration fits.




## NEXT GENERATION IMPRESSION SCANNING

Conceived with 'fully automatic impression scanning' in mind, the Identica Hybrid gained its reputation as a scanner optimized for impression scanning. For the next generation of scanners, Medit's engineering team set out to increase impression scanning capability with higher resolution cameras, advanced software algorithms and new software tools to bring you the best and most productive impression scanner ever.

## FLEXIBLE MULTI-DIE

Speed up your workflow by scanning full arch or partials with dies simultaneously. Complete your work with reduced scanning steps.



## EVOLVED CONVENIENCE OF HARDWARE & SOFTWARE FOR PRODUCTIVITY

Flexible grip provides an easy, quick mounting of single dies and full arches on the flexible multi-die. This bundled accessory is designed to "hold what you need, right where you need it". Releases easily by hand with no sticky residue left on your workpiece.

## VARIOUS STRATEGIES FOR IMPLANT CASE SCANNING

Previously, multiple implant cases could not be scanned if you only had one scan body. With Identica collLab software you can use a single scan body to capture multiple implant locations. If you have enough scan bodies you don't need to scan twice to capture the base and scan body separately. The scan body can be extracted from a single scan that contains both the base and scan body which saves you time.