

DIGITAL MATERIALS
WX
DENTAL MILLING WAX


FIT, FORM & FUNCTION

A **PERFECT** RESTORATION **BEGINS** WITH THE WAX-UP

Consistently get precise, affordable, well-fitting restorations with Jensen's digital milling wax portfolio.

ALL JENSEN DENTAL MILLING WAXES FEATURE:

- Excellent handling and performance
- Stability and strength
- Clean, ash free burn-out during casting
- Smooth castings with minimal tool marks
- Can be milled very thin without cracking or chipping
- Minimal deformation of patterns for greater accuracy



GET THE PERFECT RESTORATION EVERY TIME



DIAGNOSTIC WHITE WAX

Ideal for manufacturing diagnostic case presentations, as well as, manufacturing templates for final crowns, substructure and full contour patterns for casting and pressing.

- White color for replicating life-like final results
- Semi translucent
- Clean & simple for patient to see
- Can be carved or added to
- Available in 98x16mm discs


Please Note: Jensen offers Diagnostic White add-on wax to be used in conjunction with the Diagnostic White milling wax for adding or adjusting the milled final outcome.



SOFT BEIGE WAX

Ideal for cases requiring corrections, add-ons, and occlusal refinements of substructure and full contour patterns manufactured for casting and pressing.


- Easy to carve
- Vanilla scented
- Opaque - easy to see texture, shadows and form
- Can be carved or added to
- Available in 98x16mm discs



WAX - RED

Ideal for the manufacturing of substructure and full contour patterns for casting and pressing.

- Dark color for high contrast between wax and model
- Rigid & stiff - no warping or deformation of patterns
- Details easy to see
- Available in 98x16mm discs and 15.4mm Lava Frames



WAX - GREY

Ideal for cases requiring corrections, add-ons, and occlusal refinements of substructure and full contour patterns manufactured for casting and pressing.

- Easy to carve
- Opaque - easy to see texture, shadows and form
- Can be carved or added to
- Available in 98x16mm discs

WITH JENSEN MILLING WAXES, YOU CAN QUICKLY AND CONSISTENTLY GET THE BEST FIT, FORM AND FUNCTION WITH EVERY RESTORATION PRODUCED.