



NEWS RELEASE

DEKEMA QUALITY SINTERING FURNACES FROM JENSEN DENTAL

North Haven, Conn. – New from Jensen Dental, the DEKEMA 664 iSiC is a competitively priced sintering furnace for table top CNC milling or dependable milling center backup.

All components of the DEKEMA AUSTROMAT baSiC® and DEKEMA AUSTROMAT® µSiC® firing chambers consist of highly pure materials and provide very homogeneous temperature distribution for sintering without discoloring the sintered object. Temperature is checked and controlled by two thermocouples made of platinum/platinum rhodium.

The firing chambers are especially designed in a cylindrical shape, closed from below by a vibration-free lift. About 25 units can be placed on the DEKEMA AUSTROMAT baSiC® sintering platform, and the DEKEMA AUSTROMAT® µSiC® can sinter more than 80 units at a time



About Jensen Dental

For over 35 years, Jensen Dental has helped dental laboratories and dentists create beautiful restorations more efficiently by offering a full suite of innovative solutions backed by world-class education and support. Jensen provides quality products and services including digital CAD/CAM systems and digital materials, alloys, refining, porcelain, investment, furnaces and more.

Media Contact

Jensen Dental
Ann Pellegrini
(203) 239-2090
ap@jensendental.com
www.jensendental.com