

Lava™ Plus - High Translucency Zirconia Dyeing Liquids

VITA 3D-Master® Shade Guide Conversion Table

Lava™ Plus Dyeing Liquids	VITA 3D-Master® (Best Match)
A1	1M2
A2	2M2
A3	2R2.5
A3,5	3R2.5
A4	4M2
B1	1M1
B2	2L1.5
B3	2M3
B4	3M3

Lava™ Plus Dyeing Liquids	VITA 3D-Master® (Best Match)
C1	2L1.5
C2	3L1.5
C3	4L1.5
C4	5M2
D2	2L1.5
D3	3L1.5
D4	3L2.5
W1	0M1
W3	0M3

