

Preciso™

M205

5 Axis Dental Mill



The **Preciso M205 5 Axis Mill** consistently delivers incredibly accurate & detailed restorations by combining best in class milling technology with approved long lasting tools, driven by Preciso's proven customized CAM software.

Consistent. Accurate. Reliable.



Features

- 5 Axis Simultaneous Milling
- Up to 35° Axis Rotation
- Mills 98mm Dry Discs
- Includes powerful, versatile Preciso CAM
- 300W 60,000 RPM 4 Fold Hybrid Ceramic Bearing Spindle
- Haptic (Mechanical) Touch Point
- 16 Tool Auto Changer
- Sealed working chamber protects valuable working parts
- Integrated drawer for extra burs & (6) 98mm discs
- 165lbs, 18x21x25 in. (wxdxh)